



P Series Whole Autoclavable Micropipette



121°C sterilization



Wide range



Dispense accurately

TECHNICAL DATA

Model	Capacity Range	Test Capacity	Tip	Increment	Accuracy Error	Repeat Error
P2.5	0.1 - 2.5ul	2.5ul	10ul	0.05ul	±2.00%	±1.50%
		1.25ul			±2.50%	±2.50%
		0.25ul			±1.50%	±3.00%
P10	0.5 - 10ul	10ul	10ul	0.1ul	±0.90%	±0.80%
		5ul			±1.00%	±1.00%
		1ul			±2.00%	±1.50%
P20	2 - 20ul	20ul	200ul	0.5ul	±0.80%	±0.40%
		10ul			±1.00%	±1.00%
		2ul			±2.00%	±1.50%
P50	5 - 50ul	50ul	200ul	0.5ul	±0.60%	±0.30%
		25ul			±0.90%	±0.60%
		5ul			±1.50%	±2.00%
P100	10 - 100ul	100ul	200ul	1ul	±0.80%	±0.15%
		50ul			±1.00%	±0.40%
		10ul			±2.00%	±1.50%
P200	20 - 200ul	200ul	200ul	1ul	±0.60%	±0.15%
		100ul			±0.80%	±0.30%
		20ul			±2.00%	±1.00%
P1000	100 - 1000ul	1000ul	1000ul	5ul	±0.60%	±0.20%
		500ul			±0.70%	±0.25%
		100ul			±1.00%	±0.70%
P5000	1000 - 5000ul	5000ul	5000ul	50ul	+0.50%	±0.15%
		2500ul			±0.60%	±0.30%
		1000ul			±0.70%	±0.30%
P10000	2 - 10ml	10ml	10ml	0.1ml	±0.60%	+0.20%
		5ml			+1.20%	±0.30%
		2ml			+2.00%	±0.60%